



**Laerdal**<sup>®</sup>  
helping save lives

## Meet the 2019 AHA Feedback Requirement with Laerdal Resuscitation Training Solutions

By January 31, 2019, the American Heart Association will require the use of an instrumented directive feedback device or manikin in all AHA courses that teach the skills of adult CPR. Laerdal's portfolio of CPR training products includes feedback solutions to help you meet the requirement and deliver quality CPR training for your learners.



*Resusci Anne<sup>®</sup> QCPR and SimPad<sup>®</sup> PLUS with SkillReporter<sup>™</sup>*



*Little Anne<sup>®</sup> QCPR with SkillGuide<sup>™</sup>*

	Feedback Device	SimPad <sup>®</sup> PLUS with SkillReporter	SkillGuide <sup>™</sup>	SkillReporter (PC)	CPRmeter 2	QCPR Instructor/Learner app	SimPad <sup>®</sup> SkillReporter (Discontinued)	CPRmeter (Discontinued)
Manikin	Resusci Anne <sup>®</sup> QCPR	✓	✓	✓	✓	✓	✓	✓
	Little Anne <sup>®</sup> QCPR		✓		✓	✓		✓
	Little Anne <sup>®</sup> *				✓			✓

### KEY

- Device meets **minimum** standard as defined by AHA: instrumented, auditory and/or visual feedback on rate and depth of compression
- Device meets **optimum** standard as defined by AHA: meets the minimum standard and also includes auditory and/or visual feedback on recoil and didactic(s) for proper hand placement
- Device not compatible with manikin

*Learn more about Laerdal's feedback solutions at [Laerdal.com](http://Laerdal.com)*

\*Little Anne can be upgraded to QCPR with the Little Anne QCPR Upgrade Kit.